

Application Note

ISE - Dilutions

The following is general information regarding dilutions.

Required Information:

- What concentration is the stock solution?
- What is the volume of the required standards?
- What is the volume of the volumetric flask being used for the dilution?
- What are the volumes of the Class A pipettes?

Calculation Formula:

(Concentration A) (Volume A) = (Concentration B) (Volume B)

where:

Concentration A - the concentration of the stock solution

Volume A - the amount of stock solution needed for the dilution

Concentration B - the concentration needed

Volume B - the volume of the volumetric flask used for the dilution

Example:

Concentration of stock solution: 1000 ppm (Concentration A)

New standard: 10 ppm (Concentration B)

Volumetric flask: 100 mL (Volume B)

$(1000 \text{ ppm}) (\text{Volume A}) = (10 \text{ ppm}) (100 \text{ mL})$

Volume A = 1 mL of the stock solution

So, for this example, to make a 10 ppm standard, add 1 mL of the 1000 ppm standard to a 100 mL volumetric flask. Add DI water to mark. Mix by inversion.

Types of Dilutions:

- Serial dilution – All standards are made by diluting the previous standard. Example: The 100 ppm is made from the 1000 ppm. The 10 ppm is made from the 100 ppm.
- Non-serial dilution - Where all standards are diluted from the stock solution. Example: The 100 ppm and 10 ppm are made from the 1000 ppm.

Always use serial dilution for making standards. Non-serial dilutions can be used as check values, as this will catch dilution errors.

Recommendations:

- Use Class A glassware
- Volume A to Volume B ratio should be within two factors of ten. This will reduce the dilution error
- Dilute up to the volume of the volumetric flask
- Mix well before making another dilution with this new standard