



## *Titrant Standardization*

### **General**

As the titre of the titrant solution changes or can change, for example due to the penetration of moisture in the bottle, it is recommended that the real titre value is determined with a standard containing a known amount of water.

### **Standardization with pure water**

Pure water should be administered with a calibrated micro syringe.

The amount of water injected should be estimated so that half of the volume of the titrant burette is dispensed.

If a 5 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette, inject exactly 25 µl water (25 mg).

If a 2 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette, inject exactly 10 µl water (10 mg).

### **Programming example for 5 mg h<sub>2</sub>O/ml titrant standarization**

#### **KF reagent titre**

Calibrate

Autostart: Yes  
Cal. Every: 2 days  
Reproducibility: 0.1%

#### **Standard**

Standard ID: Yes  
Standard unit: mg  
Advised amount: 25.000 mg  
Uncertainty: 0.500 mg  
H<sub>2</sub>O in std.: 100.00%  
±H<sub>2</sub>O in std.: 0.001%

#### **Parameters**

Stirring speed: 500 rpm  
Max. burette speed: 150%/min  
Min. titr. time: 00:10 (min:sec)  
Max. titr. time: 00:02 (h:min)  
Max. volume: 10 ml

### **Standardization with HYDRANAL Standard 5.00 (Water in method)**

This calibration standard can be administered with a syringe or a pipette, its water content is: 5.00 ±0.02 mg H<sub>2</sub>O/ml at 20°C

If a 5 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette, add exactly 5 ml Standard 5.00.

If a 2 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette, add exactly 2 ml Standard 5.00.

### **Programming example for 5 mg h<sub>2</sub>O/ml titrant standarization**

#### **KF reagent titre**

Calibrate  
Autostart: Yes  
Cal. Every: 2 days  
Reproducibility: 0.1%

#### **Standard**

Standard ID: Yes  
Standard unit: ml  
Advised amount: 5.000 ml  
Uncertainty: 0.050 ml  
H<sub>2</sub>O in std.: 5.000 mg/ml  
±H<sub>2</sub>O in std.: 0.020 mg/ml

#### **Parameters**

Stirring speed: 500 rpm  
Max. burette speed: 150%/min  
Min. titr time: 00:10 (min:sec)  
Max. titr time: 00:02 (h:min)  
Max. volume: 10 ml

### **Standardization with sodium tartrate dihydrate**

Sodium tartrate dihydrate does not lose or adsorb moisture, it has a water content of 15.66% ±0.05%. This standard is a fine powder and should be administered with a weighing funnel. The exact amount is determined by difference weighing.

If a 5 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette add, exactly around 160 mg sodium tartrate dihydrate.

If a 2 mg H<sub>2</sub>O/ml titrant is used in a 10 ml burette add, exactly around 64 mg sodium tartrate dihydrate.

**Programming example for a 5 mg H<sub>2</sub>O/ml titrant  
standardization**

**KF reagent titre**

Calibrate

Autostart: Yes

Cal. Every: 2 days

Reproducibility: 0.1%

**Standard**

Standard ID: Yes

Standard unit: g

Advised amount: 0.160 g

Uncertainty: 0.001 g

H<sub>2</sub>O in std.: 15.66%

±H<sub>2</sub>O in std.: 0.050%

**Parameters**

Stirring speed: 600 rpm

Max. burette speed: 150%/min

Min. titr time: 00:30 (min:sec)

Max. titr time: 00:05 (h:min)

Max. volume: 10 ml

**Note:**

With sodium tartrate dihydrate, the solvent must be renewed after each titrant standardization.