

UltraBasic Meters



UltraBasic Benchtop


- Simple, four-key operation
- Standardization at the press of a single key.
- Auto recognition of buffers
- Large easy-to-read display
- Icons show meter status
- Meter includes AC transformer and pH electrode except for "A" models which are without electrode. Deluxe Kit includes base, electrode support arm and pH/ATC electrode.

UltraBasic Portable




- Simple, four-key operation
- Standardization at the press of a single key.
- Auto recognition of buffers
- Large easy-to-read display
- Icons show meter status
- IP-67 waterproof rating
- 9V battery operation or AC (with optional transformer)
- Meter includes 9V battery and pH electrode except for "A" models which are without electrode.

200 Series Advanced Meters



- Simultaneous display of multiple modes
- Choose model for needed measurements: pH, mV, temperature, ion, and conductivity
- Easy-to-follow menus
- Bidirectional RS232, 620 data point storage with data logging
- Multi-user features include: calibration reminder, standardization delay, measure lock and data alarm
- Meter includes electrode support arm, AC transformer and high-performance, tris-compatible pH/ATC electrode except for "A" models which are without electrode.

Electrodes



Denver Instrument offers a complete line of pH, pH/ATC, ion selective and conductivity electrodes.

Call your RepLink Representative to determine the right electrode for your application.

Quick Reference Guide to Denver Electrochemistry Modes

Measurement Modes

	pH/Temp	pH/Temp/mV	pH/Temp/ mV/Condo	pH/Temp/ mV/Ion	pH/Temp/ mV/Ion/Condo
UltraBasic Benchtop Meters	UB 5	UB 10			
200 Series Advanced Benchtop Meters		Model 215	Model 220	Model 225	Model 250
UltraBasic Portable Meters	UP 5	UP 10		UP 25	

Need more information?
Call your RepLink Representative





Precision weighing instruments since 1880

1855 Blake St. Suite 201
Denver, Colorado 80202
800.321.1135
303.431.7255
303.423.4831 FAX
www.denverinstrumentusa.com



Denver Instrument Quick Reference Guide

Electrochemistry



Competitive Cross Reference

UltraBasic Benchtop Meters

Denver			Corning			Oakton			VWR			WTW			Orion			
UB-5	02-226-211	\$395														310	13-642-551	\$556
																410Aplus	13-641-150	\$445
UB-5A	02-226-212	\$355							SB21 Kit	14002-726	\$824							
UB-10	02-226-213	\$495							SB20 Kit	14002-722	\$670					320	13-642-627	\$735
																330	13-642-552	\$817
																350	13-642-628	\$945
																350	13-642-629	\$998
																420Aplus	13-641-249	\$554
UB-10A	02-226-214	\$455	Model 430	13-641-1132	\$426				SB21	14002-728	\$773	Level 1 pH	09-500-224	\$495				
									SB20	14002-724	\$567							
UB-10 Deluxe	02-226-83	\$685	Model 430 Kit	13-641-1132	\$707	pH 510	15-500-013	\$625				Level 1 pH Kit	09-500-226	\$631				
UB-10 Deluxe Bio	02-226-91	\$725				pH 510 Bio	15-500-014	\$650										

UltraBasic Portable Meters

Denver			Corning			Oakton			VWR			WTW			Orion			
UP-5	02-226-215	\$350														210Aplus Kit	13-641-151	\$359
UP-5A	02-226-216	\$310														210Aplus	13-641-152	\$349
UP-10	02-226-217	\$425	Model 313 Kit	13-641-1008	\$455	Acorn pH 6 Kit	15-500-001	\$270	SP21 Kit	14002-704	\$695	pH315i	09-500-211	\$400		230Aplus Kit	13-641-246	\$461
			Model 314 Kit	13-641-1010	\$578	pH 10 Kit	15-500-006	\$330	SP20 Kit	14002-700	\$597	pH330i	09-500-214	\$530		250Aplus Kit	13-641-247	\$589
						pH 100 Kit	15-500-009	\$435				pH340i	09-500-217	\$670		260A Kit	13-302-36	\$790
																265A Kit	13-302-39	\$982
																610 Kit	13-641-715	\$925
UP-10A	02-226-218	\$385	Model 313	13-641-1007	\$319	Acorn pH 6	15-500-000	\$180	SP21	14002-706	\$633					260A	13-302-37	\$737
			Model 314	13-641-1009	\$445	pH 10	15-500-005	\$206	SP20	14002-702	\$515					265A	13-304-40	\$950
						pH 100	15-500-008	\$312								610	13-301-160	\$575
UP-25	02-226-219	\$695	Model 315 Kit	13-641-1012	\$932	Acorn Ion 6 Kit	15-500-041	\$374	SP30I Kit	14002-709	\$798					290Aplus Kit	13-641-248	\$742
UP-25A	02-226-220	\$655	Model 315	13-641-1011	\$584	Acorn Ion 6	15-500-040	\$313	SP30I	14002-710	\$716	pH/ion 340i	09-500-257	\$982				

200 Series Benchtop Meters

Denver			Corning			Oakton			VWR			WTW			Orion			
Model 215	02-226-105	\$1,075	Model 440 Kit	13-641-1137	\$986							Level 2 pH Kit	09-500-234	\$844		520Aplus	13-641-251	\$1,026
			Model 445 Kit	13-641-1141	\$1,270											520Aplus Kit	13-641-243	\$1,151
																535A Kit	13-641-407	\$1,395
																620 Kit	13-641-714	\$1,310
Model 215 A	02-226-153	\$920	Model 440	13-641-1136	\$718							Level 2 pH	09-500-230	\$650		535A	13-641-412	\$1,195
			Model 445	13-641-1140	\$996											620	13-301-161	\$820
Model 220	02-226-109	\$1,465				pH/Con 510	15-500-017	\$870	SR40C Kit	14002-742	\$1,570	Level 1 pH/condo Kit	09-500-294	\$1,359		550A w/electrode	13-302-89	\$1,595
																555A Kit	13-641-401	\$1,795
Model 220A	02-226-157	\$1,305							SR40C	14002-744	\$1,468	Level 1 pH/condo	09-500-292	\$1,030		550A	13-302-96	\$1,525
																555A	13-641-406	\$1,645
Model 225	02-226-113	\$1,465	Model 450 Kit	13-641-1149	\$1,569	Ion 510	15-500-019	\$710	SR30I Kit	14002-746	\$1,597	Level 2 pH/ion Kit	09-500-263	\$1,079		370	13-642-553	\$1,049
			Model 455 Kit	13-641-1149	\$2,688							Level 3 pH/ion Kit	09-500-345	\$1,289		525Aplus	13-641-342	\$1,229
																525A plus Kit	13-641-344	\$1,389
																710Aplus Kit	13-641-269	\$954
																720Aplus Kit	13-641-253	\$1,378
																920Aplus Kit	13-641-255	\$2,121
Model 225 A	02-226-161	\$1,305	Model 450	13-641-1144	\$1,296				SR30I	14002-748	\$1,494	Level 2 pH/ion	09-500-261	\$890		720Aplus	13-641-242	\$1,272
			Model 455	13-641-1148	\$2,411							Level 3 pH/ion	09-500-242	\$1,165		920Aplus	13-641-224	\$2,015
																EA940	13-641-741	\$4,081
Model 250	02-226-121	\$1,885							SR601C Kit	14002-750	\$1,983	Level 3 pH/condo Kit	09-500-312	\$1,074				
Model 250A	02-226-169	\$1,725							SR601C	14002-752	\$1,879	Level 3 pH/condo	09-500-310	\$1,585				



1.800.321.1135

www.denverinstrumentusa.com