


UltraBasic Meters

UltraBasic Benchtop




- Simple, four-key operation
- Standardization at the press of a single key.
- Auto recognition of buffers
- Large easy-to-read display
- Icons show meter status
- Meter includes AC transformer and pH electrode except for "A" models which are without electrode. Deluxe Kit includes base, electrode support arm and pH/ATC electrode.

UltraBasic Portable




- Simple, four-key operation
- Standardization at the press of a single key.
- Auto recognition of buffers
- Large easy-to-read display
- Icons show meter status
- IP-67 waterproof rating
- 9V battery operation or AC (with optional transformer)
- Meter includes 9V battery and pH electrode except for "A" models which are without electrode.

200 Series Advanced Meters



- Simultaneous display of multiple modes
- Choose model for needed measurements: pH, mV, temperature, ion, and conductivity
- Easy-to-follow menus
- Bidirectional RS232, 620 data point storage with data logging
- Multi-user features include: calibration reminder, standardization delay, measure lock and data alarm
- Meter includes electrode support arm, AC transformer and high-performance, tris-compatible pH/ATC electrode except for "A" models which are without electrode.

Electrodes



Denver Instrument offers a complete line of pH, pH/ATC, ion selective and conductivity electrodes.

Call your RepLink Representative to determine the right electrode for your application.

Quick Reference Guide to Denver Electrochemistry Modes

	Measurement Modes				
	pH/Temp	pH/Temp/mV	pH/Temp/mV/Condo	pH/Temp/mV/Ion	pH/Temp/mV/Ion/Condo
UltraBasic Benchtop Meters	UB 5	UB 10			
200 Series Advanced Benchtop Meters		Model 215	Model 220	Model 225	Model 250
UltraBasic Portable Meters	UP 5	UP 10		UP 25	

Free Electrode or Rebate!

Offer extended through 2008

Model	Rebate Amount
200 Series	\$100 or Free Electrode
UP Portable	\$50 or Free Electrode
UB Benchtop	\$50 or Free Electrode



Precision laboratory instruments since 1880



1401 17th St. Suite 750
 Denver, Colorado 80202
 800.321.1135
 303.431.7255
 303.423.4831 FAX

www.denverinstrumentUSA.com



Denver Instrument Quick Reference Guide

Electrochemistry



Competitive Cross Reference

Denver	Mettler*	Orton*	Oakton*	Nova Analytical*	VWR** (Symphony Series)
UB-5	02-226-211 \$415				
UB-5A	02-226-212 \$365				
UB-10	02-226-213 \$475				
UB-10A	02-226-214 \$425				
UB-10 w/Arm & Base	02-226-86 \$595	SevenEasy S20 01-913-806 \$670		Pinnacle 530 NC9104967 \$685	
UB-10 Deluxe Bio	02-226-83 \$765	SevenEasy S20K 01-913-807 \$820	2 Star pH 13-642-604 \$690	pH 510 Bio 15-500-013 \$646	Pinnacle 540 Kit NC9084339 \$1,377
UB-10 Deluxe	02-226-83 \$765	SevenEasy S20K 01-913-807 \$820	2 Star pH 13-642-604 \$690	pH 510 15-500-013 \$646	Pinnacle 530 Kit NC9091876 \$987
UB-10	02-226-217 \$455	SevenGo pH SG2-EK01 01-910-515 \$595	Acorn pH 6 Kit 15-500-001 \$350	Model 313 Kit NC9135754 \$636	SP70P 11388340 \$494
UP-10A	02-226-218 \$405	SevenGo pH SG2 01-910-514 \$430	3 Star pH 13-642-059 \$454	pH 110 Kit 13-620-238 \$571	Model 314 Kit NC9084633 \$808
UP-10A	02-226-218 \$405	SevenGo pH SG2 01-910-514 \$430	3 Star pH 13-642-059 \$454	Acorn pH 6 15-500-000 \$231	Model 313 NC9146889 \$446
UP-25	02-226-219 \$625			pH 110 13-620-239 \$501	Model 314 NC9063581 \$620
UP-25A	02-226-220 \$575			Acorn Ion 6 Kit 15-500-041 \$481	SP80PI 11388334 \$786
UP-25A	02-226-220 \$575	4 Star pH/ISE 13-642-062 \$860		Acorn Ion 6 15-500-040 \$410	SP80PI 11388320 \$708

UltraBASIC Portable Meters

Denver	Mettler*	Orton*	Oakton*	Nova Analytical*	VWR**
Model 215	02-226-105 \$1,255	SevenMulti S40K 01-913-809 \$1,380	3 Star pH 13-645-102 \$902	Pinnacle 545 Kit NC9096202 \$1,773	SB70P 11388368 \$630
Model 215A	02-226-153 \$1,055	SevenMulti S40 01-913-808 \$1,220	3 Star pH 13-642-052 \$567	Pinnacle 545 NC9076268 \$1,428	SB70P 11388354 \$546
Model 220	02-226-109 \$1,715	SevenMulti S30K 01-913-801 \$1,130		Pinnacle 545 15-500-017 \$943	SC70C 11388366 \$880
Model 220A	02-226-157 \$1,515	SevenEasy S30 01-913-800 \$905	3 Star Condo 13-642-054 \$782	Pinnacle 541 NC9099280 \$1,143	SB70C 11388352 \$673
Model 225	02-226-113 \$1,645	SevenMulti S70 01-913-802 \$1,660	4 Star Condo 13-645-057 \$1,210		SB80PC 11388360 \$1,015
Model 225A	02-226-161 \$1,445	SevenMulti S80K 01-913-811 \$2,080			SB80PI 11388348 \$915
Model 250	02-226-121 \$2,145	SevenMulti S80 01-913-810 \$1,920	720Aplus 13-641-242 \$1,596		
Model 250A	02-226-169 \$1,945		5 Star 13-642-058 \$1,315		SB90M5 11388342 \$1,114

200 Series Benchtop Meters

Denver	Mettler*	Orton*	Oakton*	Nova Analytical*	VWR**
Model 215	02-226-105 \$1,255	SevenMulti S40K 01-913-809 \$1,380	3 Star pH 13-645-102 \$902	Pinnacle 545 Kit NC9096202 \$1,773	SB70P 11388368 \$630
Model 215A	02-226-153 \$1,055	SevenMulti S40 01-913-808 \$1,220	3 Star pH 13-642-052 \$567	Pinnacle 545 NC9076268 \$1,428	SB70P 11388354 \$546
Model 220	02-226-109 \$1,715	SevenMulti S30K 01-913-801 \$1,130		Pinnacle 545 15-500-017 \$943	SC70C 11388366 \$880
Model 220A	02-226-157 \$1,515	SevenEasy S30 01-913-800 \$905	3 Star Condo 13-642-054 \$782	Pinnacle 541 NC9099280 \$1,143	SB70C 11388352 \$673
Model 225	02-226-113 \$1,645	SevenMulti S70 01-913-802 \$1,660	4 Star Condo 13-645-057 \$1,210		SB80PC 11388360 \$1,015
Model 225A	02-226-161 \$1,445	SevenMulti S80K 01-913-811 \$2,080			SB80PI 11388348 \$915
Model 250	02-226-121 \$2,145	SevenMulti S80 01-913-810 \$1,920	720Aplus 13-641-242 \$1,596		
Model 250A	02-226-169 \$1,945		5 Star 13-642-058 \$1,315		SB90M5 11388342 \$1,114

*Orion, Mettler, Ockton and Nova Analytical pricing via FisherSci.com, December 2007

**VWR pricing via VWR.com, December 2007



1.800.321.1135

www.denverinstrumentUSA.com

Ask About Our Free Electrode Promotion!

Offer extended through 2008