

Powder Safety Workstation From Denver Instrument

**Purchase an Aspen 200 Powder Safety Workstation
between January 1 and June 30, 2008
and receive a \$250 American Express Gift Cheque®**

Aspen 200 Features:

- Designed specifically for the containment of powder
- Saves energy by returning clean air to the room
- Constructed of polypropylene, with blue base for easy powder detection
- Includes pre-filter and HEPA filter with 99.997% efficiency at 0.3 microns
- Low-airflow and back pressure alarms warn of any potential problems
- Available in 32" or 48" widths



Please send rebate to:

Name

Company

Mailing Address

City

State

Zip

Phone

E-mail Address

To receive your rebate, please fax or mail
this form with a copy of your invoice to :

 **DENVER
INSTRUMENT**
Precision laboratory instruments since 1880

1401 17th St. Suite 750
Denver, Colorado 80202
800.321.1135
303.431.7255
303.423.4831 FAX
www.denverinstrumentUSA.com

Offer valid on purchases between January 1 and
June 30, 2008 Offer void where prohibited. Must
redeem by July 31, 2008. Allow 4 weeks for
delivery. Not valid with any other offer.