



# XI-Series High-Capacity Scales



- Auto-zero tracking
- Full digital calibration
- Full tare
- Simple, push button kg/lb conversion
- Rechargeable 6V battery or AC powered (adapter included)
- Versatile display can be pedestal or wall mounted
- Rugged, die-cast aluminum construction
- Selectable sample size in counting mode (10, 20, 50, 100)
- Adjustable display module pivots for best viewing angle
- Large, 1" (25mm) LCD numerals
- Stainless steel platform
- Low battery warning
- Overload protection
- Automatic shutoff
- RS232 interface option
- One-year warranty



## Specifications

Model	XI-30K	XI-60K	XI-120K
Readability	30kg	60kg	120kg
Division	5g	10g	20g
Platform Dimensions	14" x 18" (350mm x 450mm)		
Operating Temperature	0°C - 40°C		
Power Supply	AC Adapter, 9V/500mA		
Overall Dimensions	25.4" x 14" x 30.8" (635mm x 350mm x 770mm)		
Net Weight	27.5 lbs. (12.5kg)		
Fisher Cat. No.	02-226-200	02-226-201	02-226-202
Serial Interface Option	02-226-203		



Denver Instrument Company

*Precision laboratory instruments since 1880*

6542 Fig Street • Arvada, Colorado 80004  
 303-431-7255 • 1-800-321-1135 • Fax 303-423-4831  
[www.denverinstrument.com](http://www.denverinstrument.com)